

PDV6000plusmetal/reagent matrix

Metal	Application Note No.	Electrode set	Electrode set part No.	Standard	Electrolyte	Plate	Electrode conditioning	Other reagent	Electrolyte diluent	Cell
Ag	AN056	Bare Carbon	SP-600-000ELS-02	R-300-20STAG-01	R-300-500CLA-01	NONE	R-300-025ECS-01	R-300-025EDT-03	R-300-500ELD-01	STD
As(III)	AN029	Au film	SP-600-000ELS-05	R-300-20STAS-01	R-300-500EAS-01	R-300-025PLA-AU	R-300-025ECS-02	R-300-500EAS-02	R-300-500ELD-02	STD
As TOTAL	AN061	Au film	SP-600-000ELS-05	R-300-20STAS-01	R-300-500EAS-01	R-300-025PLA-AU	R-300-025ECS-02	NONE	R-300-500ELD-02	STD
As TOTAL	AN027	Solid Au	SP-600-000ELS-01	R-300-20STAS-01	R-300-500EAS-01	NONE	R-300-025ECS-02	NONE	R-300-500ELD-02	STD
Au	AN046	Bare Carbon	SP-600-000ELS-02	R-300-20STAU-01	R-300-500EAU-01	NONE	R-300-025ECS-01	NONE	R-300-500ELD-01	STD
Bi	AN052	Hg film	SP-600-000ELS-04	R-300-20STBI-01	R-300-500HME-01	R-300-025PLA-HG	R-300-025ECS-01	NONE	R-300-500ELD-01	STD
Cd	AN026 or AN028	Hg film	SP-600-000ELS-04	R-300-20STAA-01 or R-300-20STCD-01	R-300-500CLA-01	R-300-025PLA-HG	R-300-025ECS-01	NONE	R-300-500ELD-01	STD
Co	AN047	Cr film	SP-600-000ELS-03	R-300-20STCR-01	R-300-500CLA-01	R-300-025PLA-CR	R-300-025ECS-01	R-300-025HMA-01	R-300-500ELD-01	STD
Co & Ni	AN066	Bi film	SV LabCell	R-300-20STCO-01 Or R-300-20STNI-01	R-300-500ENI-01	R-300-500PLA-BI	NONE	R-300-500ENI-02	NONE	SV LabCell
Cr(VI)	AN032	Bare Carbon	SP-600-000ELS-02	R-300-020STCR-01	R-300-500EHG-01	NONE	R-300-025ECS-01	NONE	R-300-500ELD-01	STD
Cr(VI)	AN063	Bi film	SV LabCell	R-300-25STCR-01	R-300-500CLA-01	R-300-500PLA-BI	NONE	R-300-25CUAM-01	R-300-500ELD-01	SV LabCell
Cu	AN028	Hg film	SP-600-000ELS-04	R-300-20STAA-01 or R-300-20STCU-01	R-300-500CLA-01	R-300-025PLA-HG	R-300-025ECS-01	NONE	R-300-500ELD-01	STD
Fe	AN033	Cr film	SP-600-000ELS-03	R-300-20STFE-01	R-300-500CLA-01	R-300-025PLA-CR	R-300-025ECS-01	R-300-025HMA-01	R-300-500ELD-01	STD
Hg	AN031	Bare Carbon	SP-600-000ELS-02	R-300-20STHG-01	R-300-500EHG-01	NONE	R-300-025ECS-01	NONE	R-300-500ELD-01	STD
Hg	AN045	Au film	SP-600-000ELS-05	R-300-20STHG-01	R-300-500EHG-02	R-300-025PLA-AU	R-300-025ECS-01	NONE	R-300-500ELD-01	STD
Mn	AN048	Hg film	SP-600-000ELS-04	R-300-20STMN-01	R-300-500CLA-01	R-300-025PLA-HG	R-300-025ECS-01	NONE	R-300-500ELD-01	STD
Mo	AN059	Bi film	SV LabCell	R-300-20STMO-01	R-300-500EMO-01	R-300-500PLA-BI	NONE	NONE	NONE	SV LabCell
Ni	AN030	Cr film	SP-600-000ELS-03	R-300-20STNI-01	R-300-500CLA-01	R-300-025PLA-CR	R-300-025ECS-01	R-300-025HMA-01	R-300-500ELD-01	STD
Ni	AN060	Bi Film	SV LabCell	R-300-20STNI-01	R-300-500ENI-01	R-300-500PLA-BI	NONE	R-300-500ENI-02	R-300-500ELD-01	SV LabCell
Pd	AN051	Bare Carbon	SP-600-000ELS-02	R-300-20STPD-01	R-300-500DAR-01	NONE	R-300-025ECS-01	NONE	NONE	STD
Pb	AN026 or AN028	Hg film	SP-600-000ELS-04	R-300-20STAA-01 or R-300-20STPB-01	R-300-500CLA-01	R-300-025PLA-HG	R-300-025ECS-01	NONE	R-300-500ELD-01	STD
Sb(III)	AN053	Hg film	SP-600-000ELS-04	R-300-20STSB-01	R-300-500ESB-01	R-300-025PLA-HG	R-300-025ECS-01	R-300-500ESB-02	NONE	STD
Se	AN055	Au film	SP-600-000ELS-05	R-300-20STSE-01	R-300-500ESE-01	R-300-025PLA-AU	R-300-025ECS-01	R-300-20STCU-01	NONE	STD
Sn	AN054	Hg film	SP-600-000ELS-04	R-300-20STSN-01	R-300-500ESN-01	R-300-025PLA-HG	R-300-025ECS-01	R-300-500ESN-02	NONE	STD
Te	AN049	Au film	SP-600-000ELS-05	R-300-20STTE-01	R-300-500ESE-01	R-300-025PLA-AU	R-300-025ECS-01	NONE	NONE	STD
Tl	AN050	Hg film	SP-600-000ELS-04	R-300-20STTL-01	R-300-500CLA-01	R-300-025PLA-HG	R-300-025ECS-01	R-300-025EDT-02	R-300-500ELD-01	STD
U	AN064	Bi	SV LabCell	R-300-20STUR-01	R-300-500CLA-01	R-300-500PLA-BI	NONE	R-300-25CUAM-01	R-300-500ELD-01	SV LabCell
Zn	AN026	Hg film	SP-600-000ELS-04	R-300-20STBB-01 or R-300-20STZN-01	R-300-500CLA-01	R-300-025PLA-HG	R-300-025ECS-01	R-300-025GRA-01	R-300-500ELD-01	STD